

For Research Use Only

vwf Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66682-5-PBS



Basic Information

Catalog Number:

66682-5-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG25578

GenBank Accession Number:

GeneID (NCBI):

7450

UNIPROT ID:

P04275

Full Name:

von Willebrand factor

Purification Method:

Protein G Magarose purification

CloneNo.:

3F2E10

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51009-5, vwf Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66682-5-PBS. Detection antibody: 66682-2-PBS. Standard: Ag25578. Range: 0.098-100 ng/mL.



Cytometric bead array standard curve of MP51009-7, vwf Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66682-5-PBS. Detection antibody: 66682-6-PBS. Standard: Ag25578. Range: 0.195-100 ng/mL.