

For Research Use Only

UCP2 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66700-2-PBS



Basic Information

Catalog Number:

66700-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG19361

GenBank Accession Number:

BC011737

GeneID (NCBI):

7351

UNIPROT ID:

P55851

Full Name:

uncoupling protein 2 (mitochondrial,
proton carrier)

Calculated MW:

33 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2A1C12

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50586-1, UCP2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66700-2-PBS. Detection antibody: 66700-3-PBS. Standard: Ag19361. Range: 0.098-100 ng/mL.



Cytometric bead array standard curve of MP50586-2, UCP2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66700-2-PBS. Detection antibody: 66700-4-PBS. Standard: Ag19361. Range: 0.781-100 ng/mL.