

For Research Use Only

# CD21 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66701-2-PBS



## Basic Information

Catalog Number:

66701-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG18058

GenBank Accession Number:

BC136394

GeneID (NCBI):

1380

ENSEMBL Gene ID:

ENSG00000117322

UNIPROT ID:

P20023

Full Name:

complement component (3d/Epstein  
Barr virus) receptor 2

Calculated MW:

1092 aa, 119 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2B8A4

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP51164-1, CD21 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66701-2-PBS. Detection antibody: 66701-3-PBS. Standard: Ag18058. Range: 0.391-50 ng/mL



Cytometric bead array standard curve of MP51164-3, CD21 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66701-2-PBS. Detection antibody: 66701-6-PBS. Standard: Ag18058. Range: 0.781-100 ng/mL