

For Research Use Only

CD21 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 66701-3-PBS



Basic Information

Catalog Number: 66701-3-PBS	GenBank Accession Number: BC136394	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 1380	CloneNo.: 2G1B7
Source: Mouse	ENSEMBL Gene ID: ENSG00000117322	
Isotype: IgG2a	UNIPROT ID: P20023	
Immunogen Catalog Number: AG18058	Full Name: complement component (3d/Epstein Barr virus) receptor 2	
	Calculated MW: 1092 aa, 119 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage


Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51164-1, CD21 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66701-2-PBS. Detection antibody: 66701-3-PBS. Standard:Ag18058. Range: 0.391-50 ng/mL