

For Research Use Only

CD21 Monoclonal antibody, PBS Only (Detector)

Catalog Number:66701-6-PBS



Basic Information

Catalog Number:

66701-6-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG18058

GenBank Accession Number:

BC136394

GeneID (NCBI):

1380

ENSEMBL Gene ID:

ENSG00000117322

UNIPROT ID:

P20023

Full Name:

complement component (3d/Epstein Barr virus) receptor 2

Calculated MW:

1092 aa, 119 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

1H7A4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51164-3, CD21 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66701-2-PBS. Detection antibody: 66701-6-PBS. Standard:Ag18058. Range: 0.781-100 ng/mL.