

For Research Use Only

CHCHD4 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66718-2-PBS



Basic Information

Catalog Number:

66718-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG15339

GenBank Accession Number:

BC033775

GeneID (NCBI):

131474

UNIPROT ID:

Q8N4Q1

Full Name:

coiled-coil-helix-coiled-coil-helix domain containing 4

Calculated MW:

142 aa, 16 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

1A10A3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50491-1, CHCHD4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66718-2-PBS. Detection antibody: 66718-3-PBS. Standard: Ag15339. Range: 0.391-12.5 ng/mL