

For Research Use Only

MRPS27 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66724-2-PBS



Basic Information

Catalog Number: 66724-2-PBS	GenBank Accession Number: BC030521	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 23107	CloneNo.: 2E12E7
Source: Mouse	UNIPROT ID: Q92552	
Isotype: IgG2a	Full Name: mitochondrial ribosomal protein S27	
Immunogen Catalog Number: AG11047	Calculated MW: 168 aa, 19 kDa, 48 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Cytometric bead array standard curve of MP50532-1, MRPS27 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66724-2-PBS. Detection antibody: 66724-3-PBS. Standard:Ag11047. Range: 0.391-100 ng/mL



Cytometric bead array standard curve of MP50532-2, MRPS27 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66724-2-PBS. Detection antibody: 66724-4-PBS. Standard:Ag11047. Range: 0.781-50 ng/mL



Cytometric bead array standard curve of MP50532-3, MRPS27 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66724-2-PBS. Detection antibody: 66724-5-PBS. Standard:Ag11047. Range: 3.125-50 ng/mL