

For Research Use Only

ATG5 Monoclonal antibody, PBS Only (Capture/Detector)



Catalog Number: 66744-3-PBS

Basic Information

Catalog Number:

66744-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG16685

GenBank Accession Number:

BC002699

GeneID (NCBI):

9474

UNIPROT ID:

Q9H1Y0

Full Name:

ATG5 autophagy related 5 homolog
(S. cerevisiae)

Calculated MW:

32 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

5H8B4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50323-1, ATG5 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66744-2-PBS. Detection antibody: 66744-3-PBS. Standard: Ag16685. Range: 3.125-100 ng/mL.



Cytometric bead array standard curve of MP50323-2, ATG5 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66744-3-PBS. Detection antibody: 66744-4-PBS. Standard: Ag16685. Range: 3.125-100 ng/mL.