

For Research Use Only

Androgen receptor Monoclonal antibody, PBS Only (Detector)

Catalog Number: 66747-5-PBS



Basic Information

Catalog Number: 66747-5-PBS	GenBank Accession Number: BC132975	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 367	CloneNo.: 2E1B8
Source: Mouse	UNIPROT ID: P10275	
Isotype: IgG1	Full Name: androgen receptor	
Immunogen Catalog Number: AG17291	Calculated MW: 914 aa, 99 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51133-2, Androgen receptor Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66747-4-PBS. Detection antibody: 66747-5-PBS. Standard: Ag17291. Range: 1.563-200 ng/mL.



Cytometric bead array standard curve of MP51133-3, Androgen receptor Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66747-6-PBS. Detection antibody: 66747-5-PBS. Standard: Ag17291. Range: 0.391-200 ng/mL.