

For Research Use Only

TRAIL Monoclonal antibody, PBS Only (Detector)



Catalog Number: 66756-2-PBS

Basic Information

Catalog Number:

66756-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG25746

GenBank Accession Number:

BC032722

GeneID (NCBI):

8743

UNIPROT ID:

P50591

Full Name:

tumor necrosis factor (ligand)
superfamily, member 10

Calculated MW:

281 aa, 33 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1B10B8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50341-1, TRAIL Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66756-1-PBS. Detection antibody: 66756-2-PBS. Standard: Ag25746. Range: 3.125-200 ng/mL.