

For Research Use Only

# CD52 Monoclonal antibody, PBS Only



Catalog Number: 66784-1-PBS

## Basic Information

Catalog Number:

66784-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG17788

GenBank Accession Number:

BC000644

GeneID (NCBI):

1043

UNIPROT ID:

P31358

Full Name:

CD52 molecule

Calculated MW:

6 kDa

Observed MW:

20 kDa

Purification Method:

Protein A purification

CloneNo.:

1B10E8

## Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



human spleen tissue were subjected to SDS PAGE followed by western blot with 66784-1-Ig (CD52 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 66784-1-PBS in a different storage buffer formulation.