

For Research Use Only

CCR6 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 66801-3-PBS

Basic Information

Catalog Number:

66801-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG24195

GenBank Accession Number:

BC037960

GeneID (NCBI):

1235

UNIPROT ID:

P51684

Full Name:

chemokine (C-C motif) receptor 6

Calculated MW:

42 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

3E2A8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50358-1, CCR6 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66801-2-PBS. Detection antibody: 66801-3-PBS. Standard:Ag24195. Range: 6.25-200 ng/mL