

For Research Use Only

PRX3 Monoclonal antibody, PBS Only (Capture/Detector)



Catalog Number: 66810-2-PBS

Basic Information

Catalog Number: 66810-2-PBS	GenBank Accession Number: BC007062	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 10935	CloneNo.: 1G9A6
Source: Mouse	UNIPROT ID: P30048	
Isotype: IgG1	Full Name: peroxiredoxin 3	
Immunogen Catalog Number: AG1062	Calculated MW: 27 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50218-1, PRX3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66810-2-PBS. Detection antibody: 66810-3-PBS. Standard: Ag1062. Range: 0.781-100 ng/mL



Cytometric bead array standard curve of MP50218-2, PRX3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66810-1-PBS. Detection antibody: 66810-2-PBS. Standard: Ag1062. Range: 0.781-100 ng/mL