

For Research Use Only

PRX3 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 66810-3-PBS

Basic Information

Catalog Number:

66810-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG1062

GenBank Accession Number:

BC007062

GeneID (NCBI):

10935

UNIPROT ID:

P30048

Full Name:

peroxiredoxin 3

Calculated MW:

27 kDa

Purification Method:

Protein A purification

CloneNo.:

1F1A9

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50218-1, PRX3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66810-2-PBS. Detection antibody: 66810-3-PBS. Standard: Ag1062. Range: 0.781-100 ng/mL.