

For Research Use Only

PRDX1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66820-4-PBS



Basic Information

Catalog Number:

66820-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG8821

GenBank Accession Number:

BC007063

GeneID (NCBI):

5052

UNIPROT ID:

Q06830

Full Name:

peroxiredoxin 1

Calculated MW:

199 aa, 22 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1B11C4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50445-2, PRDX1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66820-4-PBS. Detection antibody: 66820-5-PBS. Standard: Ag8821. Range: 0.781-100 ng/mL.