

For Research Use Only

IBA1 Monoclonal antibody, PBS Only (Capture/Detector)

Catalog Number: 66827-5-PBS



Basic Information

Catalog Number: 66827-5-PBS	GenBank Accession Number: BC009474	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 199	CloneNo.: 1D7A7
Source: Mouse	ENSEMBL Gene ID: ENSG00000204472	
Isotype: IgG1	UNIPROT ID: P55008	
Immunogen Catalog Number: AG28236	Full Name: allograft inflammatory factor 1	
	Calculated MW: 17 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51023-2, IBA1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66827-3-PBS. Detection antibody: 66827-5-PBS. Standard:Ag28236. Range: 1.563-100 ng/mL.



Cytometric bead array standard curve of MP51023-3, IBA1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66827-5-PBS. Detection antibody: 66827-4-PBS. Standard:Ag28236. Range: 3.125-100 ng/mL.