

For Research Use Only

IL-12RB1/CD212 Monoclonal antibody, PBS Only (Capture)

Catalog Number:66829-2-PBS



Basic Information

Catalog Number: 66829-2-PBS	GenBank Accession Number: BC029121	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 3594	CloneNo.: 1F11F3
Source: Mouse	UNIPROT ID: P42701	
Isotype: IgG1	Full Name: interleukin 12 receptor, beta 1	
Immunogen Catalog Number: AG6920	Calculated MW: 662 aa, 73 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50569-1, IL-12RB1/CD212 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66829-2-PBS. Detection antibody: 66829-3-PBS. Standard:Ag6920. Range: 0.098-100 ng/mL



Cytometric bead array standard curve of MP50569-2, IL-12RB1/CD212 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66829-2-PBS. Detection antibody: 66829-4-PBS. Standard:Ag6920. Range: 0.098-50 ng/mL