

For Research Use Only

ASGR2 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66831-2-PBS



Basic Information

Catalog Number: 66831-2-PBS	GenBank Accession Number: BC017251	Purification Method: Protein A Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 433	CloneNo.: 1G7F12
Source: Mouse	UNIPROT ID: P07307	
Isotype: IgG2b	Full Name: asialoglycoprotein receptor 2	
Immunogen Catalog Number: AG2037	Calculated MW: 39 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50757-1, ASGR2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66831-2-PBS. Detection antibody: 66831-3-PBS. Standard:Ag2037. Range: 0.781-25 ng/mL



Cytometric bead array standard curve of MP50757-2, ASGR2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66831-2-PBS. Detection antibody: 66831-4-PBS. Standard:Ag2037. Range: 0.781-25 ng/mL