

For Research Use Only

# AGER/RAGE Monoclonal antibody, PBS Only (Detector)

Catalog Number: 66833-4-PBS



## Basic Information

Catalog Number:

66833-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG9657

GenBank Accession Number:

BC020669

GeneID (NCBI):

177

UNIPROT ID:

Q15109

Full Name:

advanced glycosylation end product-specific receptor

Calculated MW:

404 aa, 43 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

3A7B4

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP51094-2, AGER/RAGE Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66833-2-PBS. Detection antibody: 66833-4-PBS. Standard:Ag9657. Range: 0.391-12.5 ng/mL.