

For Research Use Only

G6PC Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66860-2-PBS



Basic Information

Catalog Number:

66860-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG17839

GenBank Accession Number:

BC130478

GeneID (NCBI):

2538

UNIPROT ID:

P35575

Full Name:

glucose-6-phosphatase, catalytic subunit

Calculated MW:

357 aa, 40 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

3D7E6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50826-1, G6PC Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66860-2-PBS. Detection antibody: 66860-3-PBS. Standard: Ag17839. Range: 0.098-100 ng/mL.