

For Research Use Only

G6PC Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66860-4-PBS



Basic Information

Catalog Number:
66860-4-PBS

Size:
1 mg/ml

Source:
Mouse

Isotype:
IgG2a

Immunogen Catalog Number:
AG17839

GenBank Accession Number:
BC130478

GeneID (NCBI):
2538

UNIPROT ID:
P35575

Full Name:
glucose-6-phosphatase, catalytic
subunit

Calculated MW:
357 aa, 40 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:
2F5H3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50826-2, G6PC Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66860-4-PBS. Detection antibody: 66860-1-PBS. Standard: Ag17839. Range: 0.098-100 ng/mL.