

For Research Use Only

# G6PC Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66860-5-PBS



## Basic Information

Catalog Number:

66860-5-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG17839

GenBank Accession Number:

BC130478

GeneID (NCBI):

2538

UNIPROT ID:

P35575

Full Name:

glucose-6-phosphatase, catalytic  
subunit

Calculated MW:

357 aa, 40 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

3C12D3

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50826-3, G6PC Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66860-5-PBS. Detection antibody: 66860-6-PBS. Standard: Ag17839. Range: 0.098-6.25 ng/mL