

For Research Use Only

CD8A Monoclonal antibody, PBS Only (Capture)



Catalog Number: 66868-2-PBS

Basic Information

Catalog Number: 66868-2-PBS	GenBank Accession Number: BC025715	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 925	CloneNo.: 1E12B3
Source: Mouse	ENSEMBL Gene ID: ENSG00000153563	
Isotype: IgG1	UNIPROT ID: P01732	
Immunogen Catalog Number: AG11111	Full Name: CD8a molecule	
	Calculated MW: 235 aa, 26 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50246-1, CD8A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66868-2-PBS. Detection antibody: 66868-3-PBS. Standard: Ag11111. Range: 3.125-100 ng/mL.