

For Research Use Only

# CD8A Monoclonal antibody, PBS Only (Detector)



Catalog Number: 66868-3-PBS

## Basic Information

Catalog Number:

66868-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG11111

GenBank Accession Number:

BC025715

GeneID (NCBI):

925

ENSEMBL Gene ID:

ENSG00000153563

UNIPROT ID:

P01732

Full Name:

CD8a molecule

Calculated MW:

235 aa, 26 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1B2A4

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50246-1, CD8A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66868-2-PBS. Detection antibody: 66868-3-PBS. Standard:Ag11111. Range: 3.125-100 ng/mL.