

For Research Use Only

STIL Monoclonal antibody, PBS Only (Capture)



Catalog Number: 66876-2-PBS

Basic Information

Catalog Number: 66876-2-PBS	GenBank Accession Number: BC126223	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 6491	CloneNo.: 1G4A3
Source: Mouse	UNIPROT ID: Q15468	
Isotype: IgG1	Full Name: SCL/TAL1 interrupting locus	
Immunogen Catalog Number: AG19124	Calculated MW: 1287 aa, 143 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50266-1, STIL Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66876-2-PBS. Detection antibody: 66876-3-PBS. Standard:Ag19124. Range: 1.563-200 ng/mL



Cytometric bead array standard curve of MP50266-2, STIL Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66876-2-PBS. Detection antibody: 66876-4-PBS. Standard:Ag19124. Range: 1.563-200 ng/mL