

For Research Use Only

STIL Monoclonal antibody, PBS Only (Detector)



Catalog Number: 66876-3-PBS

Basic Information

Catalog Number:

66876-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG19124

GenBank Accession Number:

BC126223

GeneID (NCBI):

6491

UNIPROT ID:

Q15468

Full Name:

SCL/TAL1 interrupting locus

Calculated MW:

1287 aa, 143 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

1B4F1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50266-1, STII Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66876-2-PBS. Detection antibody: 66876-3-PBS. Standard: Ag19124. Range: 1.563-200 ng/mL.