

For Research Use Only

TRIM45 Monoclonal antibody

Catalog Number: 66896-1-Ig



Basic Information

Catalog Number:

66896-1-Ig

Size:

1700 µg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG28300

GenBank Accession Number:

BC034943

GeneID (NCBI):

80263

UNIPROT ID:

Q9H8W5

Full Name:

tripartite motif-containing 45

Calculated MW:

562aa,62 kDa; 580aa,64 kDa

Observed MW:

62-64 kDa

Purification Method:

Protein G purification

CloneNo.:

3A7H11

Recommended Dilutions:

WB 1:1000-1:6000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human

Positive Controls:

WB : HeLa cells, HepG2 cells, Jurkat cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 66896-1-Ig (TRIM45 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.