

For Research Use Only

Integrin alpha 6 Monoclonal antibody, PBS Only



Catalog Number: 66906-1-PBS

Basic Information

Catalog Number: 66906-1-PBS	GenBank Accession Number: BC136456	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 3655	CloneNo.: 2B12A1
Source: Mouse	UNIPROT ID: P23229	
Isotype: IgG2b	Full Name: integrin, alpha 6	
Immunogen Catalog Number: AG25866	Calculated MW: 1073 aa, 119 kDa	
	Observed MW: 115-120 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



A431 cells were subjected to SDS PAGE followed by western blot with 66906-1-Ig (Integrin alpha-6 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 66906-1-PBS in a different storage buffer formulation.