

For Research Use Only

# Integrin alpha-6 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 66906-3-PBS



## Basic Information

Catalog Number:

66906-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG25866

GenBank Accession Number:

BC136456

GeneID (NCBI):

3655

UNIPROT ID:

P23229

Full Name:

integrin, alpha 6

Calculated MW:

1073 aa, 119 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1A2C2

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50396-1, Integrin alpha-6 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66906-2-PBS. Detection antibody: 66906-3-PBS. Standard: Ag25866. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP50396-3, Integrin alpha-6 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66906-5-PBS. Detection antibody: 66906-3-PBS. Standard: Ag25866. Range: 1.563-100 ng/mL.