

For Research Use Only

Integrin alpha-6 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66906-4-PBS



Basic Information

Catalog Number: 66906-4-PBS	GenBank Accession Number: BC136456	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 3655	CloneNo.: 2F3B3
Source: Mouse	UNIPROT ID: P23229	
Isotype: IgG1	Full Name: integrin, alpha 6	
Immunogen Catalog Number: AG25866	Calculated MW: 1073 aa, 119 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50396-2, Integrin alpha-6 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66906-4-PBS. Detection antibody: 66906-2-PBS. Standard: Ag25866. Range: 3.125-100 ng/mL.