

For Research Use Only

MAVS; VISA Monoclonal antibody, PBS Only (Capture/Detector)

Catalog Number:66911-3-PBS



Basic Information

Catalog Number: 66911-3-PBS	GenBank Accession Number: BC044952	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 57506	CloneNo.: 1B4G6
Source: Mouse	UNIPROT ID: Q7Z434	
Isotype: IgG1	Full Name: mitochondrial antiviral signaling protein	
Immunogen Catalog Number: AG5949	Calculated MW: 57 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50509-1, MAVS; VISA Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66911-2-PBS. Detection antibody: 66911-3-PBS. Standard: Ag5949. Range: 0.391-100 ng/mL



Cytometric bead array standard curve of MP50509-2, MAVS; VISA Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66911-3-PBS. Detection antibody: 66911-4-PBS. Standard: Ag5949. Range: 0.098-100 ng/mL