

For Research Use Only

# MAVS; VISA Monoclonal antibody, PBS Only (Detector)

Catalog Number:66911-4-PBS



## Basic Information

Catalog Number:

66911-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG5949

GenBank Accession Number:

BC044952

GeneID (NCBI):

57506

UNIPROT ID:

Q7Z434

Full Name:

mitochondrial antiviral signaling

protein

Calculated MW:

57 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2A9B7

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50509-2, MAVS; VISA Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66911-3-PBS. Detection antibody: 66911-4-PBS. Standard: Ag5949. Range: 0.098-100 ng/mL



Cytometric bead array standard curve of MP50509-3, MAVS; VISA Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66911-1-PBS. Detection antibody: 66911-4-PBS. Standard: Ag5949. Range: 0.391-50 ng/mL