

For Research Use Only

CD155/PVR Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66913-2-PBS



Basic Information

Catalog Number:

66913-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG28303

GenBank Accession Number:

BC015542

GeneID (NCBI):

5817

UNIPROT ID:

P15151

Full Name:

poliovirus receptor

Calculated MW:

45 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

1G8B8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50402-1, CD155/PVR Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66913-2-PBS. Detection antibody: 66913-3-PBS. Standard: Ag28303. Range: 3.125-100 ng/mL.



Cytometric bead array standard curve of MP50402-2, CD155/PVR Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66913-2-PBS. Detection antibody: 66913-4-PBS. Standard: Ag28303. Range: 3.125-100 ng/mL.