

For Research Use Only

TFG Monoclonal antibody, PBS Only (Detector)

Catalog Number: 66916-2-PBS



Basic Information

Catalog Number: 66916-2-PBS	GenBank Accession Number: BC023599	Purification Method: Euglobulin precipitation
Size: 1 mg/ml	GeneID (NCBI): 10342	CloneNo.: 4D1D6
Source: Mouse	UNIPROT ID: Q92734	
Isotype: IgM	Full Name: TRK-fused gene	
Immunogen Catalog Number: AG27697	Calculated MW: 400 aa, 43 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50894-1, TFG Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66916-1-PBS. Detection antibody: 66916-2-PBS. Standard: Ag27697. Range: 0.391-200 ng/mL.