

For Research Use Only

BST2 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66919-2-PBS



Basic Information

Catalog Number: 66919-2-PBS	GenBank Accession Number: BC033873	Purification Method: Protein A Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 684	CloneNo.: 3C6D10
Source: Mouse	UNIPROT ID: Q10589	
Isotype: IgG2a	Full Name: bone marrow stromal cell antigen 2	
Immunogen Catalog Number: AG4430	Calculated MW: 180 aa, 20 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50437-1, BST2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66919-2-PBS. Detection antibody: 66919-3-PBS. Standard: Ag4430. Range: 0.391-100 ng/mL.



Cytometric bead array standard curve of MP50437-2, BST2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66919-2-PBS. Detection antibody: 66919-4-PBS. Standard: Ag4430. Range: 1.563-100 ng/mL.