

For Research Use Only

Integrin beta-4 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66922-3-PBS



Basic Information

| | | |
|---------------------------------------------|----------------------------------------------|----------------------------------------------------------------|
| Catalog Number: 66922-3-PBS | GenBank Accession Number: BC126411 | Purification Method: Protein G Magarose purification |
| Size: 1 mg/ml | GeneID (NCBI): 3691 | CloneNo.: 3A12B4 |
| Source: Mouse | UNIPROT ID: P16144 | |
| Isotype: IgG1 | Full Name: integrin, beta 4 | |
| Immunogen Catalog Number: AG17395 | Calculated MW: 1752 aa, 195 kDa | |

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51169-1, Integrin beta-4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66922-3-PBS. Detection antibody: 66922-4-PBS. Standard: Ag17395. Range: 0.391-100 ng/mL