

For Research Use Only

CFTR Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66928-2-PBS



Basic Information

Catalog Number:

66928-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG27810

GenBank Accession Number:

NM_000492

GeneID (NCBI):

1080

UNIPROT ID:

P13569

Full Name:

cystic fibrosis transmembrane
conductance regulator (ATP-binding
cassette sub-family C, member 7)

Calculated MW:

168 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1D2B10

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50929-1, CFTR Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66928-2-PBS. Detection antibody: 66928-3-PBS. Standard: Ag27810. Range: 0.391-25 ng/mL