

For Research Use Only

NTF3 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66933-2-PBS



Basic Information

Catalog Number: 66933-2-PBS	GenBank Accession Number: BC069773	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 4908	CloneNo.: 2B10D9
Source: Mouse	UNIPROT ID: P20783	
Isotype: IgG1	Full Name: neurotrophin 3	
Immunogen Catalog Number: AG12790	Calculated MW: 257 aa, 29 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50760-1, NTF3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66933-2-PBS. Detection antibody: 60520-1-PBS. Standard: Ag12790. Range: 0.195-200 ng/mL.