

For Research Use Only

MRRF Monoclonal antibody, PBS Only (Detector)

Catalog Number: 66942-3-PBS



Basic Information

Catalog Number:

66942-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG28583

GenBank Accession Number:

BC013049

GeneID (NCBI):

92399

UNIPROT ID:

Q96E11

Full Name:

mitochondrial ribosome recycling
factor

Calculated MW:

262 aa, 29 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1C11F3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50707-1, MRRF Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66942-2-PBS. Detection antibody: 66942-3-PBS. Standard:Ag28583. Range: 0.781-50 ng/mL



Cytometric bead array standard curve of MP50707-2, MRRF Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66942-4-PBS. Detection antibody: 66942-3-PBS. Standard:Ag28583. Range: 1.563-100 ng/mL