

For Research Use Only

ART4 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 66949-1-PBS

Basic Information

Catalog Number: 66949-1-PBS	GenBank Accession Number: BC074727	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 420	CloneNo.: 2A11E3
Source: Mouse	UNIPROT ID: Q93070	
Isotype: IgG1	Full Name: ADP-ribosyltransferase 4 (Dombrock blood group)	
Immunogen Catalog Number: AG28173	Calculated MW: 314 aa, 36 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50324-2, ART4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66949-1-PBS. Detection antibody: 66949-4-PBS. Standard: Ag28173. Range: 1.563-100 ng/mL.