

For Research Use Only

# ITGB3 Monoclonal antibody, PBS Only (Capture/Detector)

Catalog Number: 66952-2-PBS



## Basic Information

Catalog Number:

66952-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG27834

GenBank Accession Number:

BC127666

GeneID (NCBI):

3690

UNIPROT ID:

P05106

Full Name:

integrin, beta 3 (platelet glycoprotein

IIIa, antigen CD61)

Calculated MW:

788 aa, 87 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

3C9G3

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50447-1, ITGB3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66952-2-PBS. Detection antibody: 66952-3-PBS. Standard: Ag27834. Range: 0.391-100 ng/mL



Cytometric bead array standard curve of MP50447-2, ITGB3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66952-4-PBS. Detection antibody: 66952-2-PBS. Standard: Ag27834. Range: 3.125-100 ng/mL