

For Research Use Only

ITGB3 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 66952-3-PBS



Basic Information

Catalog Number: 66952-3-PBS	GenBank Accession Number: BC127666	Purification Method: Protein A Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 3690	CloneNo.: 3D5A2
Source: Mouse	UNIPROT ID: P05106	
Isotype: IgG2a	Full Name: integrin, beta 3 (platelet glycoprotein IIIa, antigen CD61)	
Immunogen Catalog Number: AG27834	Calculated MW: 788 aa, 87 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage


Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50447-1, ITGB3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66952-2-PBS. Detection antibody: 66952-3-PBS. Standard: Ag27834. Range: 0.391-100 ng/mL.