

For Research Use Only

NOTUM Monoclonal antibody, PBS Only (Detector)

Catalog Number: 66956-3-PBS



Basic Information

Catalog Number: 66956-3-PBS	GenBank Accession Number: BC060882	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 147111	CloneNo.: 1B7H5
Source: Mouse	UNIPROT ID: Q6P988	
Isotype: IgG2a	Full Name: notum pectinacylesterase homolog (Drosophila)	
Immunogen Catalog Number: AG6338	Calculated MW: 56 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50968-1, NOTUM Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66956-2-PBS. Detection antibody: 66956-3-PBS. Standard:Ag6338. Range: 0.098-100 ng/mL.