

For Research Use Only

LXN Monoclonal antibody, PBS Only (Detector)



Catalog Number: 66957-2-PBS

Basic Information

Catalog Number: 66957-2-PBS	GenBank Accession Number: BC005346	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 56925	CloneNo.: 1G2A4
Source: Mouse	UNIPROT ID: Q9BS40	
Isotype: IgG2a	Full Name: latexin	
Immunogen Catalog Number: AG27922	Calculated MW: 222 aa, 26 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP50023-1, LXN
Monoclonal Matched Antibody Pair, PBS Only.
Capture antibody: 66957-1-PBS. Detection
antibody: HRP-conjugated 66957-2-PBS. Standard:
Ag27922. Range: 0.313-20 ng/mL.