

For Research Use Only

# HK2 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66974-4-PBS



## Basic Information

Catalog Number: 66974-4-PBS	GenBank Accession Number: BC021116	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 3099	CloneNo.: 1H2B7
Source: Mouse	UNIPROT ID: P52789	
Isotype: IgG1	Full Name: hexokinase 2	
Immunogen Catalog Number: AG16895	Calculated MW: 102 kDa	

## Applications

Tested Applications:  
Cytometric bead array, Indirect ELISA

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -80°C.  
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50448-2, HK2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66974-4-PBS. Detection antibody: 66974-5-PBS. Standard: Ag16895. Range: 1.563-100 ng/mL.



Cytometric bead array standard curve of MP50448-3, HK2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66974-4-PBS. Detection antibody: 66974-1-PBS. Standard: Ag16895. Range: 0.781-100 ng/mL.