

For Research Use Only

NFKB1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66992-2-PBS



Basic Information

Catalog Number: 66992-2-PBS	GenBank Accession Number: BC051765	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 4790	CloneNo.: 1A1F1
Source: Mouse	UNIPROT ID: P19838	
Isotype: IgG2b	Full Name: nuclear factor of kappa light polypeptide gene enhancer in B-cells 1	
Immunogen Catalog Number: AG5832	Calculated MW: 105 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51119-1, NFKB1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66992-2-PBS. Detection antibody: 66992-1-PBS. Standard: Ag5832. Range: 3.125-100 ng/mL.