

For Research Use Only

ICAM4 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67014-3-PBS



Basic Information

Catalog Number:

67014-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG27628

GenBank Accession Number:

BC029364

GeneID (NCBI):

3386

UNIPROT ID:

Q14773

Full Name:

intercellular adhesion molecule 4
(Landsteiner-Wiener blood group)

Calculated MW:

271 aa, 29 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

1F11D7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50429-1, ICAM4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67014-2-PBS. Detection antibody: 67014-3-PBS. Standard: Ag27628. Range: 0.098-100 ng/mL.