

For Research Use Only

# STBD1 Monoclonal antibody, PBS Only (Capture/Detector)



Catalog Number: 67018-2-PBS

## Basic Information

Catalog Number:

67018-2-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG28613

GenBank Accession Number:

BC022301

GeneID (NCBI):

8987

UNIPROT ID:

O95210

Full Name:

starch binding domain 1

Calculated MW:

41 kDa

Purification Method:

Protein G purification

CloneNo.:

1B3B3

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50130-1, STBD1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67018-2-PBS. Detection antibody: 67018-3-PBS. Standard:Ag28613. Range: 1.563-100 ng/mL



Cytometric bead array standard curve of MP50130-2, STBD1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67018-1-PBS. Detection antibody: 67018-2-PBS. Standard:Ag28613. Range: 1.563-100 ng/mL