

For Research Use Only

NFS1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67021-4-PBS



Basic Information

Catalog Number:

67021-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG8017

GenBank Accession Number:

BC018471

GeneID (NCBI):

9054

UNIPROT ID:

Q9Y697

Full Name:

NFS1 nitrogen fixation 1 homolog (S. cerevisiae)

Calculated MW:

50 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

1A9E1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50834-2, NFS1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67021-4-PBS. Detection antibody: 67021-5-PBS. Standard:Ag8017. Range: 0.781-200 ng/mL.