

For Research Use Only

# NFS1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67021-5-PBS



## Basic Information

Catalog Number:

67021-5-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG8017

GenBank Accession Number:

BC018471

GeneID (NCBI):

9054

UNIPROT ID:

Q9Y697

Full Name:

NFS1 nitrogen fixation 1 homolog (S.

cerevisiae)

Calculated MW:

50 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1F6G9

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50834-2, NFS1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67021-4-PBS. Detection antibody: 67021-5-PBS. Standard:Ag8017. Range: 0.781-200 ng/mL.